



## THE EFFECTS

1. Reduces the pain
2. Decrease the inflammatory process
3. Absorbs edema
4. Accelerates tissue scarring
5. Accelerates the recovery of activity




*A real professional tool for a personal use*



## USE TIPS

3 sessions / days following surgery

2 sessions / days for the treatment of pain and inflammation



easyCryo<sup>®</sup>  
MOOVE Rent



Taking charge of the return and home care of patients.



[www.easycryo.fr](http://www.easycryo.fr)



Zone ECOPARC  
27 Rue Mezière Christin  
34690 FABRÈGUES - France

Tél. +33 (0)4 99 64 21 05  
Fax. +33 (0)4 99 64 21 06  
olivia.beldame@easycryo.fr






## THE ADVANTAGES

- + 20 to 30 minutes care
- + Works with water chilled and / or ice cubes
- + Use on mains and / or on battery (5 hours)
- + Machine washable splints protective covers



## THE SPLINTS




ANKLE



ELBOW / WRIST



SHOULDER



HIP



THIGH



KNEE



## TECHNICAL DETAILS

### ICECOMPRESSION

	MOOVE RENT
Number of exits	1
Min / Max compression	0 mmHg / 75 mmHg
Compression	alternating / constant
Water temperature	1° à 8°C recommended
Program	1 (alterable)
Cycle of Cold and Rest	Function repeat
CE med	LNE/GMED 0459



## PATIENT MANAGEMENT

The Moove Rent model offers compressive cryotherapy, a recommended care technique in orthopedic surgery. This therapy combines ice water cryotherapy and air compression to reduce pain by two, accelerates healing and offers a faster recovery of physical activity.

In addition to an analgesic treatment, this treatment method will promote a better recovery of patients

