



## THE COMPRESSIVE CRYOTHERAPY

IceCompression is an innovative device combining Cryotherapy and Compression at the same time

The compressive cryotherapy is thus strongly recommended in orthopedics, sports and rheumatology- rehabilitation in the following indications: inflammation - trauma - acute and recurrent pain - edema - drainage - tendinopathies - recovery.

**Recommendations: 1 to 3 times a day for 30 minutes. For better recovery, it is advisable to maintain the effects of cold as soon as you leave the clinic.**

This cryotherapy meets the requirements of outpatient surgery and is an essential adjunct to the optimization of care.

## THE EFFECTS

**REDUCES** the pain

**BLOCKS** the inflammatory process

**ABSORBS** edemas

**ACCELERATES** tissue scarring

Treatment 1/1		
0:30:00	Exit	
Water Automati...	Splint Mask	Air Dynamic 1/30 cycles 30s/30s
0°C/8°C		75 mmHg
6°C		

## THE CONCEPT OF ICECOMPRESSION

This medical device makes it possible to treat the daily pains, the traumas and the consequences of an operation in orthopedic surgery.

Refrigerated water is sent into a splint (specific for each area of the body). The refrigerated water can be combined with maintained or alternating pressure, which will accentuate all the benefits of the cold and will limit the medical drugs intake.

A care protocol is pre-set into the software and all settings are adjustable depending on the patient needs.

The mobile water tank can also be placed in the freezer to increase the cryo capacity. The pre-established program allows you to start the treatment very quickly.



## THE ADVANTAGES

- Programmable and modifiable system
- Chilled water tank: with ice cubes or at the freezer (3h cold guarantee)
- Automatic purging of splints for a better hygiene.
- Machine washable splints protective covers
- Compatible with the dental mask after tooth, jaw or face surgery.



**THE SPLINTS**



ELBOW / WRIST



MASK



THIGH (L. 50CM)

**UPPER LIMB**

- Shoulder right and left
- Elbow / wrist

**FACE**

- Mask

**LOWER LIMB**

- Knee (L. 38 cm)
- Thigh (L. 50 cm)
- Ankle
- Hip right and left

**TECHNICAL DETAILS**

	MOOVE
Number of exits	1
Min / max compression	0 mmHg / 75 mmHg
Compression	dynamic/constant
Water temperature	1° à 8°C recommended
Program	1 (alterable)
Cycle of Cold and Rest	Function repeat
CE med	LNE/GMED 0459

HEALTH PROFESSIONALS  
AND HOME CARE PROVIDERS



The Moove



Zone ECOPARC  
27 Rue Mezière Christin  
34690 FABRÈGUES - France

Tél. +33 (0)4 99 64 21 05  
Fax. +33 (0)4 99 64 21 06  
olivia.beldame@easycryo.fr



www.easycryo.fr

MKS Paris : Sarl au capital de 10 000€ - SIRET : 51822767300028 / N° Intracom : FR 11518227673

